

Blood Bag Tube Sealer



An ISO 9001:2015 & CE Certified Company

Arnav scientific equipment bench top sealer ATS-02 is designed to meet high standards in blood safety. ATS-02 helps to keep the blood bag system completely sterile. It ensures quality sealing with clear rupture line for easy tear-apart segments. Consistent performance, safety features & monitoring system makes it ideally suited for easy & long lasting usage.

SALIENT FEATURES:

- Senses the blood bag tube for automatic sealing process
- Anti-spark and overheat protection system to avoid damage to system and tubing
- Visual indication for:-
 - Ready to seal
 - Seal process
 - Power supply
- Easy to remove sealing head cap for quick cleaning
- Seal with clear break line for snap-apart separation
- Splash guard for technician safety.



+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Blood Bag Tube Sealer



An ISO 9001:2015 & CE Certified Company

Particulars	ATS-01
Model No.	ATS-01
Dimension (W x D x H) (mm)	175X225X125 mm
Weight (Kg)	04
RF Frequency	40.68 MHZ
RF Output	80 watt
Seal time	1.5 Sec (auto adjustable)
Sealing Type	Automatic
Input Voltage	115/230 V, 50/60 Hz

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com