

Blood Bag Tube Stripper



An ISO 9001:2015 & CE Certified Company

The blood bag tube stripper is a “plier-like” tool, intended to be used with PVC blood tubing, to strip undiluted blood from the donor tubing. The blood bag tube stripper is easy to use and highly ergonomic.



SALIENT FEATURES:

- This tube stripper can be used with different sizes of blood bag tubes due to the screw mechanism that allows to adjust the distance between the rollers.
- Its lightweight body makes it extremely comfortable for everyday use.
- The exterior structure composition makes the blood bag tube stripper resistant to scratches and rust.
- **Material:** Aluminum
- **Surface Finish:** Chrome



+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com