


# Blood Bag Weight Balance




An ISO 9001:2015 & CE Certified Company

Arnav Blood Bag Weight Balance Is designed to provide accurate weight of blood bag carrier/bucket prior to centrifugation in order to avoid imbalance centrifugation operation.

## SALIENT FEATURES:



- The actual weight the contents is displayed in real time and can be adjusted simultaneously adding on removing weights.
- The differential weight in grams is also displayed at. The same time Double weight scale increases the productivity in blood banks by reducing the time required during balancing components separation process in bank centrifuge.
- Metal Body Powder coated.
- LCD display.
- Auto calibration.
- Stainless steel top pan.
- External power supply.



+91 8080007145 / 9973706825

## About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,  
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,  
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ [arnavscientificequipment@gmail.com](mailto:arnavscientificequipment@gmail.com)  
[sales@arnavscientificequipment.com](mailto:sales@arnavscientificequipment.com)

🌐 [www.arnavscientificequipment.com](http://www.arnavscientificequipment.com)  
[www.arnavscientific.com](http://www.arnavscientific.com)

# Blood Bag Weight Balance



An ISO 9001:2015 & CE Certified Company

Particulars	AWS-02
Dimension in mm (W x D x H)	300 x 350 x 150
Weight	5 Kg
Display	LCD 16x2
Feet	Adjustable
Accuracy	+/- 1gm
Transducer	Load Cell
Stainless Steel Pan Dimension	130x220
Power Supply	230 AC Volts 50Hz.

## Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,  
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,  
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ [arnavscientificequipment@gmail.com](mailto:arnavscientificequipment@gmail.com)  
[sales@arnavscientificequipment.com](mailto:sales@arnavscientificequipment.com)

🌐 [www.arnavscientificequipment.com](http://www.arnavscientificequipment.com)  
[www.arnavscientific.com](http://www.arnavscientific.com)