

# Blood Bank Refrigerator



An ISO 9001:2015 & CE Certified Company

Blood bank refrigerator (BBR) is a designated temperature controlled refrigeration machine specially designed to store blood bags. As the name indicates, it is widely used in blood banks and hospitals.

## SALIENT FEATURES:

ASE micro controller based Vertical Blood Bank Refrigerator with digital temperature Data, inner chamber made of stainless steel (304), outer body of mild steel duly powder coated, temp up to 2 to 6°C, with built in Micro controlled based temperature indicator cum controller having audio visual alarm. having low and High Temperature audio visual alarm. And Door open Alarm, Power failure Audio Visual alarm Facility available - accuracy of +/-0.5°C. Refrigerant HFC CFC Free. Power supply single phase 220volt 50Hz.



+91 8080007145 / 9973706825

## About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,  
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,  
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ [arnavscientificequipment@gmail.com](mailto:arnavscientificequipment@gmail.com)  
[sales@arnavscientificequipment.com](mailto:sales@arnavscientificequipment.com)

🌐 [www.arnavscientificequipment.com](http://www.arnavscientificequipment.com)  
[www.arnavscientific.com](http://www.arnavscientific.com)

# Blood Bank Refrigerator



An ISO 9001:2015 & CE Certified Company

Particulars	ABR-100	ABR-200	ABR-300	ABR-400
Temperature	2°C to 6°C	2°C to 6°C	2°C to 6°C	2°C to 6°C
Inner Size (L, D, H)	450X450X500	450X500X900	550X500X1100	550X560X1300
External Size (L, D, H)	570X650X920	620X620X1320	670X700X1550	720X780X1670
Storage Volume	100Ltrs	200Ltrs	300Ltrs	400Ltrs
Storage Bags	60 Bags	120 Bags	180 Bags	240 Bags
No. of Basket	03	04	05	06
Tray & Chamber MOC	SS-304	SS-304	SS-304	SS-304

## Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,  
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,  
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ [arnavscientificequipment@gmail.com](mailto:arnavscientificequipment@gmail.com)  
[sales@arnavscientificequipment.com](mailto:sales@arnavscientificequipment.com)

🌐 [www.arnavscientificequipment.com](http://www.arnavscientificequipment.com)  
[www.arnavscientific.com](http://www.arnavscientific.com)