

Blood Donor Couch



An ISO 9001:2015 & CE Certified Company

Blood Donor Couch is designed to provide a comfortable position for the donor. All surface of the couch are smooth, easily accessible and can be kept clean.

SALIENT FEATURES:



The advantage of a blood donor couch by providing a one piece lying surface with more flexibility in its function. Besides your patients requirements for a comfortable sitting position it also provides a very relaxing lying position.

Facility for blood collection from both sides. Provides a Comfortable position for the donor. Variable positioning for either arm with comfortable wide arm rests. Tilt adjustment can be done by using remote control. Comfortable Chair type with soft padding for cushioning and rexin cover. Two geared motor. Ensuring safety and comfort to the donor 4 Heavy Duty Lockable castor wheels for easy mobility.

+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Blood Donor Couch



An ISO 9001:2015 & CE Certified Company

Particulars	ADC-02
Power Supply	230 VAC
Weight	60 KG
Lifting Capacity	200 KG
Control	Corded Remote
Power Consumption	100W
Movement Actuation	DC Motor
Length of Head and Back Rest	80cm
Chair Total Width From Inside	60cm
Length of the Seat	60cm
Length of the Feet Rest	50cm
Height of the chair base on the floor	55cm

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com