

Plasma Expressor



An ISO 9001:2015 & CE Certified Company

Plasma expressor is a benchtop unit designed to extract blood components from the blood bags. It is a manual system which doesn't require any electricity and due to its strong structure, least maintenance is required. It is widely used in blood banks, hospitals and research institutes.

SALIENT FEATURES:

- Manual operations, no electricity required.
- Strong design to withstand pressure.
- Stainless steel body to prevent corrosion.
- Leveling feet at bottom for convenient usage.
- Rubber knob on handle for proper handling.
- Complies with CE, ISO standards.
- Easy to use.
- Transparent acrylic plates used to observe the separation.
- Different volumes of blood bag can be used for separation.

APPLICATIONS:

- It is used to separate blood components in hospitals, laboratories, research institutes and blood banks.



+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Plasma Expressor



An ISO 9001:2015 & CE Certified Company

Particulars	APE - 01
Model No.	APE - 01
Clamping Force	20 N
External Dimension (W x D x H)	165x250x250 mm
Separates	Plasma, Platelet
Applicable Blood Bag	100 ml, 200 ml, 300 ml, 400 ml
Package Dimension (W x D x H)	310x310x240 mm
Net Weight	2.5 Kg
Gross Weight	2.8 Kg

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com