

-40 Plasma Freezer



Arnav scientific equipment manufactures Plasma Freezer that provides an ideal freezing environment for the preservation of Plasma, and test samples.

An ISO 9001:2015 & CE Certified Company



SALIENT FEATURES:

ASE micro controller based Vertical Plasma freezer with digital temperature Data, inner chamber made of stainless steel (304), outer body of mild steel duly powder coated, temp up to -40°C , with built in Micro controlled based temperature indicator cum controller having audio visual alarm having low and High Temperature audio visual alarm. And Door open Alarm, Power failure Audio Visual alarm Facility available -40°C in AC room having an accuracy of $\pm 3^{\circ}\text{C}$. Refrigerant HFC CFC Free. Power supply single phase 220volt. 50Hz.

| Model No. | APF-120 | APF-200 | APF-300 |
|--------------------------------|--------------|--------------|--------------|
| Inner Size (H x D x W) (mm) | 430x540x510 | 550x550x660 | 650x650x720 |
| External Size (H x D x W) (mm) | 1270x630x690 | 1390x640x840 | 1490x740x900 |
| No. of Basket | 2 | 3 | 4 |

+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

+91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com