

Plasma Thawing Bath



An ISO 9001:2015 & CE Certified Company

Arnav Scientific Equipment thawing bath offer quick and optimum thawing of plasma at required temperature levels. Available in two variants Refrigerated & Non-Refrigerated.


SALIENT FEATURES:

APTB-10, is a refrigerated thawing bath for preparation of cryoprecipitate at required temperature of 37°C

APTB-10, is a non-refrigerated thawing bath for quick thawing of plasma at required temperature of 37°C.



- Microprocessor based precise temperature control.
- High performance 1.5 kw heater & compressor to achieve the set temp.
- Over temperature safety limit, audio visual alarm.
- LED monitor to display the working status of thawing bath.
- Continuous water circulation for temp uniformity & rapid thawing.
- Audio visual alarm after elapsed time.
- Capacity to thaw 15 plasma bags.



+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Plasma Thawing Bath



An ISO 9001:2015 & CE Certified Company

Particulars	APTB-10	
Temperature	+37°C (accuracy +/-0.2°C)	+ 37°C (accuracy +/-0.2°C)
External dimension w x d x h (mm)	910 x 600 x 1090	760 x 530 x 640
Thawing timer	Adjustable up to 2 Hrs	
Internal body material	Stainless steel 304 Grade	
External body material	Powder coated CRCA Steel	
Noise level	Less than 65 db (A)	
Recommended line voltage corrector	3 KVA	
Electrical	220-240 volts,50 Hz, single phase	

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com