

Platelet Incubator



An ISO 9001:2015 & CE Certified Company

Micro controller based platelet incubators are designed to store Platelet concentrates in continues motion at controlled temperature.

SALIENT FEATURES:



ASE micro controller based Vertical PLATELET INCUBATOR with digital temperature Data, inner chamber made of stainless steel (304), outer body of mild steel duly powder coated, temp up to 10 to 22°C, with built in Micro controlled based temperature indicator cum controller having audio visual alarm. having low and High Temperature audio visual alarm. And Door open Alarm, Power failure Audio Visual alarm Facility available - accuracy of +/-0.5°C Refrigerant HFC CFC Free. Power supply single phase 220volt 50Hz.

+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Platelet Incubator



An ISO 9001:2015 & CE Certified Company

Particulars	API-125
Temperature	10°C to 22°C
Inner Size (L x D x H)	600 x 500 x 450
External Size (L x D x H)	570 x 650 x 920
Storage Volume	125 Ltrs
Storage Blood Bags	48 Bags
No. of Basket	08
Tray & Chamber MOC	SS-304

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com