

Refrigerated Centrifuge




An ISO 9001:2015 & CE Certified Company

Micro controller based laboratory centrifuge is used in pathology laboratories and research centres. Motorised Safety door interlock on the centrifuge ensures that the rotor does not run with door open while the cover cannot be opened unless rotor comes to complete stop.

SALIENT FEATURES:

- Program Status (1 - 10)
- Cooling On/ Off Status
- Parameters
- Set Value Status
- Process Value Status
- Acceleration Status
- Deceleration Status
- RCF Value
- Current Status (door Open/ In Process / Process Complete/ Process Terminated/ Imbalance)
- Start Key
- Stop Key
- Program Key
- Increment / Decrement Key
- Enter Key
- Fast Cool Key (4°C & 22°C)
- Door Open Key



Processing Display



6 Bag



8 Bag



12 Bag

+91 8080007145 / 9973706825

About Us

Arnav Scientific Equipment is a Manufacturer of the Blood Bank Processing Equipment, Blood Bank Storage Equipment, Blood Bank Collection Equipment. We are manufacturing these products using premium grade raw material that is procured from the authentic vendors of the market. We offer these products at reasonable rates and deliver these within the promised time-frame.

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate,
Behind Hyundai Showroom, Waliv Phata, Sativali Rd.,
Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com

Refrigerated Centrifuge



An ISO 9001:2015 & CE Certified Company

TECHNICAL SPECIFICATION

Particulars	ARC - 06	ARC - 08	ARC - 12
Width	758 mm	758 mm	758 mm
Depth	990 mm	990 mm	990 mm
Height	1015 mm	1015 mm	1015 mm
Input Power	220-240V AC	220-240V AC	220-240V AC
Temp. Range	4°C to 22°C	4°C to 22°C	4°C to 22°C
Weight	370 Kg approx.	400 Kg approx.	430 Kg approx.
Recommended	5 KVA	7.5 KVA	7.5 KVA

ROTOR SPECIFICATION

Model No.	Rotor Details			
ARC - 06	6 Place aluminum swing out rotor head to hold 6 metal carrier having capacity to hold 1 bags with plastic bucket & wind shield 4 mm thick	4400	6000	1000 ml x 6
ARC - 08	4 Place Aluminum swing out rotor head to hold 4 metal carrier having capacity to hold 2 Bags with Plastic bucket & wind shield of 4 mm thick	4400	6125	1000 ml x 8
ARC - 12	6 Place Aluminum swing out rotor head to hold 6 metal carrier having capacity to hold 2 Bags with Plastic bucket wind shield of 4 mm thick	4400	6385	1000 ml x 12

Our Product Range

- Blood Bank Storage Equipment
- Blood Bank Processing Equipment
- Blood Bank Collection Equipment
- Laboratory Equipment
- Blood Bank Equipment
- Blood Bank Refrigerator
- Blood Donor Chair
- Laboratory Refrigerator
- Refrigerated Blood Bank Centrifuge
- Low Temperature Freezer
- Ice Lined Refrigerator
- Stability Chamber

📍 Corp Off.: Gala No. 19, Lusi Industrial Estate, Behind Hyundai Showroom, Waliv Phata, Sativali Rd., Vasai (E), Palghar 401208, Maharashtra, India.

☎ +91 8080007145 / 9973706825 / 7743838033

✉ arnavscientificequipment@gmail.com
sales@arnavscientificequipment.com

🌐 www.arnavscientificequipment.com
www.arnavscientific.com